NOTICE D'INSTRUCTION
POUR LES VIS AUTO-COMPRESSIFS STÉRILES DYNASTYSTEM®
INSTRUCTION LEAFLET FOR DYNASTYSTEM®
STERILE SELF-COMPRESSIVE SCREWS

IFU-DYN-001 Version 3
Du 2022-04-26

Fabriqué en France par
Made in France by
NEOSTEO
Mallève 2A
1, boulevard Jean Moulin
44100 NANTES – France
Tel. : +33 (0)2 36 56 96 70
Fax. : +33 (0)2 51 70 61 34
E-mail : contact@neosteo.com
Website: www.neosteo.com

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English Instructions
INSTRUCTION LEAFLET FOR NEOSTEO DYNASTYSTEM® STERILE SELF-COMPRESSIVE SCREWS & SNAP-OFF SELF-COMPRESSIVE SCREWS

	Manufacturer: NeosteO SAS Mallève 2A 1 boulevard Jean Moulin 44100 NANTES France	Tel. : +33 (0)2 36 56 96 70 Fax : +33 (0)2 51 70 61 34 E-mail : contact@neosteo.com Website: www.neosteo.com SIREN : 514 138 270	
Product reference	Batch number: XXX-YY	YYMM-YYMM Expiry date: YYMM-YYMM	YYMM-YYMM Manufacture date: YYMM-YYMM

Caution: Federal law restricts this device to the sale, distribution and the use by or on the order of a physician.

NEOSTEO will not accept any responsibility for failure to observe all of the following instructions.

1. General Information
Before using the devices, the surgeon must read and understand the contents of this product leaflet and the information specific to each device. This information is available from the distributor or manufacturer (operating techniques).
The devices can only be implanted in public or private hospital establishments, in a department of orthopedic surgery, and only by surgeons who have been trained and have sufficient experience in the field of osteotomy.
This device is for single use only. Re-use of an explanted or modified device is prohibited.

2. Description
An implant from the DynastSystem® range of self-compressive screw consists of a self-compressive screw.
The complete list of product names and references in the range are available in the dedicated surgical technique.
All of the parts of the implant are made of Ti6Al4V-EL titanium alloy for implants (standard ISO 5832-3 and ASTM F136 for Ti6Al4V).

Fixation principle
The osteosynthesis screw allows the compression of two bone fragments. The screw pitch difference between the head thread and the screw end thread compresses the two fragments, the monobloc screw ensures the mechanical resistance of the assembly. The important height of the distal thread ensures an excellent anchoring into the spongy bone, thus improving the implant performances and the osteosynthesis efficiency.

IMPORTANT : In order to achieve best performance from the system always seek to obtain optimal hold of the screws in bone, particularly using the near and far cortex in the metaphysis-diaphyseal area.

Performance of the device
Implants provide stability and anchorage on the two fixation fragments. Implants therefore enable bone to be reconstructed, providing mechanical stability for the bone whilst it is consolidating.

3. Indications
The Self-Compressive Screws are intended for the fixation of bone fractures and for bone reconstruction in the hand and in foot/surgery.
The Snap-Off Self-compressive Screws are intended for the fixation of bone fractures and for bone reconstruction in forefoot surgery.

4. Contraindications
• Acute or chronic local or systemic infections
• Absence of musculo-cutaneous coverage, severe vascular defects affecting the parts concerned
• Bone damage inconsistent with correct insertion of the screws in the bone
• Motor or neurological deficiency or behavioral disorders which risk exposing the internal fixation to abnormal mechanical forces.
• Allergy to titanium

5. Handling the product
Some implants delivered sterile are sterilized by ethylene oxide method. Implants are packaged sterile in a double bag. Ensure that the product is unpackaged under conditions and following a procedure which maintains its sterility. The implant packaging must never be damaged and the implant cannot be unpackaged or repackaged before its final use. If the sterility of any implant from the range has been lost, particularly if the package has deteriorated or if the implant comes into contact with a non-sterile surface which does not form part of the procedure it must be resterilised or used. Contact with such surfaces, the organic tissues and liquids may be a severe risk of patient contamination.
Under no circumstances will NeosteO repackage devices that have already been used.
Do not implant the device if the expiry date has passed.

6. Traceability
All of the information ensuring the traceability of the products is available in their packaging (follow up labels for patient records). This enables the implant to be monitored and assures the traceability of the product which is essential for patient safety.

7. Restrictions, ancillary equipment and related operative techniques
A NeosteO self-compressive screw and snap-off self-compressive screw cannot be used in association with other devices. Combination with other devices of different brands is strictly prohibited. The implant must not be deformed or modified in any way. The implants should only be inserted using the NeosteO ancillaries specific to them. Their full list and description are available in the surgical technique dedicated to the implants.
The different operating techniques used to position the range implants are available from product distributors or manufacturer. These should be read and understood before the implants are first used and must be followed during any surgical procedure.
The accuracy of the measuring instruments are detailed in the dedicated surgical technique.

8. Precautions
It is essential that the following points are observed during the surgical procedure.
• Ensure that the product is correct for its intended use.
• DynastSystem® screws are designed for manual insertion thanks to their double grip, which speeds up screwing. In all cases, manual insertion must be performed manually using the screwdriver provided until the head is flush with the bone surface.
• Recommended placement: **Monocortical placement only**
Canulated screws: Mono- or bicortical placement
The protruding tip can affect soft tissue if it exceeds the second cortex. In all cases, select the screw length so as to maximise bone grip, while limiting protrusion of the tip.
• Do not bend the implant during the contrast screwing.
• The wire gauge should only be used with 80 mm long and Ø1.0 mm Kirschner wires.
• To be efficient, the fixation system needs the implant to be screwed until the head is flush with the surface of the bone.
• Do not forget to complete the patient's file with the identification labels of the sterile products implanted.

9. Risks of implantation / Adverse effects
Surgical procedure of resection may involve risks of nosocomial infection, nerve compression and rejection of the implant which can be avoided by the skills and experience of the surgeon who must have been trained to perform the procedure, and by the hospital staff. The sterility of the implant must be maintained by the procedures used by the hospital staff.
The most typical and commonest effects when internal fixation implants are positioned are: infection, delayed wound healing, residual pain, under or over-correction, joint stiffness, loss of correction, nonunion, delayed bone union, traumatic and/or vascular complications due to surgery, bone fracture, hardware migration, deformation or breakage, hardware discomfort.

MRSA Safety Information:
Non-clinical testing was performed to evaluate the NeosteO Self-compressive screw. Non-clinical testing shows that the NeosteO Self-compressive screws are MR Conditional. A patient with one of these implants can be scanned safely in an MR system under the following conditions:
• Static magnetic field of 3-Tesla or less
• Minimum spatial gradient magnetic field of 1900-gauss/cm (19-T/m)
• Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode
In non-clinical testing, a NeosteO Self-compressive screw is expected to produce a maximum temperature rise of less than or equal to 6.1°C using an MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg for 15-minutes (per pulse sequence) of continuous scanning in a 1.5-Tesla MR system. The presence of NeosteO Self-compressive screws implant may produce an image artifact.

For any further information about the use of these products from the range please contact your supplier.

10. Factors liable to compromise the success of implantation
• Severe vascular defect, loss of bone viability
• Severe osteoporosis, loss of bone substance
• Very severe deformity or trauma with loss of bone substance or soft tissue
• Local bone tumor
• Systemic metabolic or genetic disorders
• Infectious diseases
• Drug addiction and/or tendency to drug and medicines abuse.
• Obesity
• Intense physical activity (e.g.: competitive sports or heavy physical work)

IMPORTANT : It is important that the patient be informed about the influence which these complications may have on the success of the procedure. It is recommended that patients be given all useful advice about measures to reduce the effects of these contraindications.

11. Post-operative protocol – Removal of the device
The post-operative procedure is the conventional protocol intended for the same surgical indications (AO protocol for osteotomies; 3 weeks minimum of partial discharge with orthopedic footwear, then normal weight-bearing of the patient).

It is recommended that the implanted material be removed as soon as no longer demonstrated to be of use (bone consolidation after the internal fixation). It is the responsibility of the surgeon to decide whether to explant the device for a smooth post-operative course following the surgical procedure.

The device can be removed using the appropriate screwdriver (consult the surgical technique manual). Disposal of the contaminated device must be carried out according to the instructions in the following paragraph.
12. Reuse – Disposal and treatment of waste
The device is single use. It must not be re-used. Amongst the usual risks likely to happen if the device is re-used :
• device contamination that can lead to infections
• oxidation or corrosion of the device that can lead to inflammation, allergies or adverse tissue effects
• deterioration /early breakage of the device that can lead to physical trauma, loss of correction or non-union.

NeosteO cannot in any situation guarantee the stability, absence of contamination and mechanical performances of the device which has already been used. All waste and residual materials must be treated as contaminated waste and removed in accordance with current regional and national legislation.

When disposing of sharp objects, it is imperative to wear suitable protective equipment for handling such objects in order to avoid any risk of superficial or deep injury, and cross-contamination.

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